Instructions For Use

SAS-MX Alkaline Hb-10

Cat. No. 100800

SAS-MX Hb-10 alkaline

Fiche technique

Réf. 100800

SAS-MX Alkalisches Hb-10

Anleitung

Kat. Nr. 100800

Alkaline Hb-10 SAS-MX

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Cod. 100800

Hb-10 Alcalina SAS-MX

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Contents

English	1
Français	7
Deutsch	13
Italiano	20
Español	26

INTENDED PURPOSE

The SAS-MX Alkaline Hb-10 kit is intended for the separation of human haemoglobins by agarose gel electrophoresis.

Haemoglobins (Hb) are a group of proteins whose chief functions are to transport oxygen from the lungs to the tissues and carbon dioxide in the reverse direction. They are composed of polypeptide chains called globin, and iron protoporphyrin haem groups. A specific sequence of amino acids constitutes each of the four polypeptide chains. Each normal haemoglobin molecule contains one pair of alpha and one pair of non-alpha chains. In normal adult haemoglobin (HbA), the non-alpha chains are called beta. The non-alpha chains of foetal haemoglobin are called gamma. A minor (3%) haemoglobin fraction called HbA₂ contains alpha and delta chains. Two other chains are formed in the embryo.

The major haemoglobin in the erythrocytes of the normal adult is HbA and there are small amounts of HbA₂ and HbF. In addition, over 400 mutant haemoglobins are now known, some of which may cause serious clinical effects, especially in the homozygous state or in combination with another abnormal haemoglobin. Wintrobe¹ divides the abnormalities of haemoglobin synthesis into three groups.

- (1) Production of an abnormal protein molecule (e.g. sickle cell anaemia).
- (2) Reduction in the amount of normal protein synthesis (e.g. thalassaemia).
- (3) Developmental anomalies (e.g. hereditary persistence of foetal haemoglobin (HPFH)).

The two mutant haemoglobins most commonly seen are HbS and HbC. Hb Lepore, HbE, HbG-Philadelphia, HbD-Los Angeles, and HbO-Arab may be seen less frequently².

Electrophoresis is generally considered the best method for separating and identifying haemoglobinopathies. The protocol for haemoglobin electrophoresis involves step-wise use of two systems^{3,8}. Initial electrophoresis is performed in alkaline buffers. However, because of the electrophoretic similarity of many structurally different haemoglobins, the evaluation must be supplemented by acid buffer electrophoresis which measures a property other than electrical charge. This method is based on the complex interactions of the haemoglobin with an alkaline electrophoretic buffer and the agarose support. The SAS-MX Alkaline Hb-10 procedure is a simple procedure requiring minute quantities of haemolysate to provide complementary evidence (along with the results from SAS-MX Acid Hb-10 analysis) of the presence of HbS, HbC and HbF as well as several other abnormal hemoglobins.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. SAS-MX Alkaline Hb Gel

Contains agarose in a Tris / EDTA / Glycine buffer with sodium azide as preservative. The gel is ready for use as packaged.

2. Tris / Borate Buffer Concentrate

Contains concentrated Tris / Borate buffer with sodium azide and thiomersal as preservatives. Dilute the contents of the bottle to 1 litre with purified water and mix well.

3. Acid Blue Stain Concentrate

Contains concentrated Acid Blue stain. Dilute the contents of the bottle to 700ml with purified water. Stir overnight and filter before use. Store in a tightly stoppered bottle.

4. Haemoglobin Lysing Agent

Contains Triton X-100 in purified water with potassium cyanide and thiomersal as preservative. The Lysing Agent is ready for use as packaged.

5. Destain Solution Concentrate

Contains concentrated Destain Solution. Dilute the contents of the bottle to 2 litres with purified water. Store in a tightly stoppered bottle.

6. Other Kit Components

Each kit contains Instructions For Use and sufficient Sample Application Templates and Blotters A and C to complete 10 gels.

STORAGE AND SHELF-LIFE

1. SAS-MX Alkaline Hb Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. Tris / Borate Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C.

3. Acid Blue Stain

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration.

4. Haemoglobin Lysing Agent

The Lysing Agent should be stored at 15...30°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

5. Destain Solution

The destain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted destain solution is stable for 6 months at 15...30°C. Cloudiness may indicate deterioration of the destain solution.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 SAS-MX Chamber

Cat. No. 1525 EPS600 Power Supply

Cat. No. 5328 AA₂ Hemo Control

Cat. No. 5329 ASA₂ Hemo Control

Cat. No. 5330 AFSA₂ Hemo Control

Cat. No. 5331 AFSC Hemo Control

Drying oven with forced air capable of 60...70°C

Fixative solution: Mix 500ml of methanol and 500ml of purified water. Add 100ml of glacial acetic acid.

Mix well and store in a tightly stoppered bottle.

Saline solution (0.85% NaCl)

SAMPLE COLLECTION AND PREPARATION

Freshly collected EDTA or heparin anticoagulated blood is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 1 week. For optimal results, saline washed red cells should be used to prepare lysates. This removes possible interference from plasma proteins.

- a) Mix 200µL of well mixed whole blood with 1000µL of saline solution.
- b) Centrifuge to sediment the red cells.
- c) Remove 1000µL of the supernatant and discard.
- d) Add a further 1000µL of saline solution and mix well.
- e) Repeat steps b-d x2.
- f) Following the final centrifugation, remove 1000µL of supernatant and treat the remaining sample as whole blood, or remove all of the supernatant and treat the remaining sample as washed packed cells.

Dilute each patient sample / control to a haemoglobin concentration of 1.0-2.0 g/dL with Haemoglobin Lysing Agent.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Blot the gel surface with a blotter C, discard the blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in step 5.
3. Apply 3µl of sample to each slit and allow to absorb for 5 minutes.
4. Whilst samples are absorbing, pour 30ml of buffer into each inner section of the SAS-MX Chamber.
5. Following sample absorption, blot the template with the blotter A retained from step 2 and remove both blotter and template.
6. Position the gel in the chamber agarose side down, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 150 volts, 30 minutes
8. At the end of electrophoresis, dry the gel at 60...70°C. **NOTE:** The drying of the gel should take no more than 5 minutes to prevent diffusion of bands. If this cannot be achieved, fix the gel in fixative solution for 5 minutes prior to drying.
9. Immerse the dry gel in stain solution for 10 minutes.

10. Destain the dry gel in 2 x 1 minute washes of destain solution or until the background is clear.
11. Dry the gel.

INTERPRETATION OF RESULTS

Qualitative Evaluation:

The possible identity of the haemoglobin types present in the samples can be determined by visual evaluation of the completed gel. The Hemo controls provide a marker for band identification.

Quantitative Evaluation:

The relative percent of each haemoglobin type on the gel can be determined by densitometry of the completed gel at 595nm.

Figure 1 shows the possible identity of the most commonly encountered haemoglobin types.


FIG. 1

Most haemoglobin variants cause no discernible clinical symptoms, so are of interest primarily to research scientists. Variants are clinically important when their presence leads to sickling disorders, thalassaemia syndromes, life long cyanosis, haemolytic anaemias or erythrocytosis, or if the heterozygote is of sufficient prevalence to warrant genetic counselling. The combinations of HbS-S, HbS-D-Los Angeles, and HbS-O Arab lead to serious sickling disorders². Several variants including HbH, E-Fort Worth and Lepore cause a thalassaemic blood picture².

The two variant hemoglobins of greatest importance in terms of frequency and pathology are HbS and HbC¹. Sickle cell anaemia (HbSS) is a cruel and lethal disease. It first manifests itself at about 5-6 months of age. The clinical course presents agonising episodes of pain and temperature elevations with anaemia, listlessness, lethargy, and infarct in virtually all organs of the body. The individual with homozygous HbCC suffers mild haemolytic anaemia which is attributed to the precipitation or crystallization of HbC within the erythrocytes. Cases of HbSC disease are characterised by haemolytic anaemia that is milder than sickle-cell anaemia.

The thalassaemias are a group of haemoglobin disorders characterised by hypochromia and microcytosis due to the diminished synthesis of one globin chain (the α or β) while synthesis of the other chain proceeds normally^{9,10}. This unbalanced synthesis results in unstable globin chains. These precipitate within the red cell, forming inclusion bodies that shorten the life span of the cell.

In α -thalassaemias, the α chains are diminished or absent, and in the β -thalassaemia, the β chains are affected. Another quantitative disorder of haemoglobin synthesis, hereditary persistent foetal haemoglobin (HPFH), represents a genetic failure of the mechanisms that turn off gamma chain synthesis at about four months after birth, which results in a continued high percentage of HbF. It is a more benign condition than the true thalassaemias and persons homozygous for HPFH have normal development, are asymptomatic and have no anemia¹⁰.

The most common haemoglobin abnormalities:

Sickle Cell Trait

This is a heterozygous state showing HbA and HbS and a normal amount of HbA₂. Results on citrate agar show haemoglobins in the HbA and HbS migratory positions (zones).

Sickle Cell Anaemia

This is a homozygous state showing almost exclusively HbS, although a small amount of HbF may also be present.

Sickle-C Disease

This is a heterozygous state demonstrating HbS and HbC.

Sickle Cell-Thalassaemia Disease

This condition shows HbA, HbF, HbS, and HbA₂.

In Sickle Cell β_0^- -Thalassaemia HbA is absent.

In Sickle Cell β_+^- -Thalassaemia HbA is present in reduced quantities.

Thalassaemia-C Disease

This condition shows HbA, HbF and HbC.

C Disease

This is a homozygous state showing almost exclusively HbC.

Thalassaemia Major

This condition shows HbF, HbA and HbA₂.

LIMITATIONS

Some abnormal haemoglobins have similar electrophoretic mobilities and must be differentiated by other methodologies.

Further testing required:

1. Acid Hb electrophoresis may be a necessary follow up test for confirmation of abnormal haemoglobins detected.
2. Globin chain analysis (both acid and alkaline) and structural studies may be necessary in order to positively identify some of the more rare haemoglobins.
3. Anion exchange column chromatography is the most accurate method for quantitating HbA₂. Helena BioSciences Sickle-Thal Quik Column Method (Cat. No. 5334) for quantitation of HbA₂ in the presence of HbS, or the Helena BioSciences Beta-Thal HbA₂ Quik Column Procedure (Cat. No. 5341) are recommended. HbA₂ quantitation is one of the most important diagnostic tests in the diagnosis of β -thalassaemia trait.
4. Low levels of HbF (1-10%) may be accurately quantitated by radial immunodiffusion using the Helena BioSciences HbF-QuikPlate Procedure (Cat. No. 9325).

REFERENCE VALUES

At birth, the majority of haemoglobin in the erythrocytes of the normal individual is foetal haemoglobin, HbF. Some of the major adult haemoglobin, HbA, and a small amount of HbA₂, are also present. At the end of the first year of life and through adulthood, the major haemoglobin present is HbA with up to 3.7% HbA₂ and less than 2% HbF.

PERFORMANCE CHARACTERISTICS

Reproducibility

	Within Gel (n=10)		Between Gel (n=100)	
	Mean (%)	CV (%)	Mean (%)	CV (%)
HbA	45.9	1.6	45.9	1.8
HbF	37.8	1.3	38.1	1.7
HbS	14.5	3.2	14.2	4.4
HbA ₂	1.8	7.8	1.9	17.8

Linearity

The linearity of the method is a function of densitometer specification as well as gel performance. It is recommended that each customer determine the linearity of the method based upon the densitometer in use in the laboratory.

BIBLIOGRAPHY

1. Wintrobe, Maxwell M., Clinical Hematology, 6th Edition, Lea and Febiger, Philadelphia, 1967. pages 145-167.
2. Fairbanks, V.F., 'The Nomenclature and Taxonomy of Hemoglobin Variants', Diagnostic Medicine, Nov/Dec., 53-58, 1980.
3. Schneider, R.G., Hightower, B.J. and Barwick, R.C., 'Laboratory Identification of the Hemoglobins', Lab Management, 1981; August: 29-43.
4. Center for Disease Control, Laboratory Methods for Detecting Hemoglobinopathies, U.S. Department of Health and Human Services/Public Health Service, 1984.
5. Schneider, R.G., 'Methods for Detection of Hemoglobin Variants and Hemoglobinopathies in the Routine Clinical Laboratory', CRC Critical Reviews in Clinical Laboratory Sciences, 1978.
6. Schneider, R.G., Hightower, B., Hosty, T.S., Ryder, H., Tomlin, G., Atkins, R., Brimhall, B., and Jones, R.T., 'Abnormal Hemoglobins in a Quarter Million People', Blood, 1976; 48(5) : 629-637.
7. Huisman, T.H.J. and Schroeder, W.A., 'New Aspects of the Structure, Function and Synthesis of Hemoglobins'. CRC Press, Cleveland, 1971.
8. Schmidt, R.M., Huisman, T.H.J., and Lehmann, H., The Detection of Hemoglobinopathies. CRC Press, Cleveland, 1974.
9. Weatherall, D.J. and Clegg, J.B., The Thalassemia Syndromes, Blackwell Scientific Publications, Oxford, 1972.
10. Lehman, H. and Huntsman, R.G., Man's Haemoglobins, J.B. Lippincott Co., Philadelphia, 1974.

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LITERATUR

- Wintrobe, Maxwell M., Clinical Hematology, 6th Edition, Lea and Febiger, Philadelphia, 1967. Seite 145-167.
- Fairbanks, V.F., 'The Nomenclature and Taxonomy of Hemoglobin Variants', Diagnostic Medicine, Nov/Dec., 53-58, 1980.
- Schneider, R.G., Hightower, B.J. and Barwick, R.C., 'Laboratory Identification of the Hemoglobins', Lab Management, 1981; August: 29-43.
- Center for Disease Control, Laboratory Methods for Detecting Hemoglobinopathies, U.S. Department of Health and Human Services/Public Health Service, 1984.
- Schneider, R.G., 'Methods for Detection of Hemoglobin Variants and Hemoglobinopathies in the Routine Clinical Laboratory', CRC Critical Reviews in Clinical Laboratory Sciences, 1978.
- Schneider, R.G., Hightower, B., Hosty, T.S., Ryder, H., Tomlin, G., Atkins, R., Brimhall, B., and Jones, R.T., 'Abnormal Hemoglobins in a Quarter Million People', Blood, 1976; 48(5) : 629-637.
- Huisman, T.H.J. and Schroeder, W.A., 'New Aspects of the Structure, Function and Synthesis of Hemoglobins'. CRC Press, Cleveland, 1971.

8. Schmidt, R.M., Huisman, T.H.J.,and Lehmann, H., The Detection of Hemoglobinopathies. CRC Press, Cleveland, 1974.
9. Weatherall, D.J. and Clegg, J.B., The Thalassemia Syndromes, Blackwell Scientific Publications, Oxford, 1972.
10. Lehman, H. and Huntsman, R.G., Man's Haemoglobins, J.B. Lippincott Co., Philadelphia, 1974.

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