## Instructions For Use

---

### SAS-I Serum Protein Cat No. 200100

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**INTENDED PURPOSE**

The SAS-I Serum Protein kit is intended for the separation and quantitation of serum proteins by agarose gel electrophoresis.

Serum contains over 100 individual proteins, each with a specific set of functions which are subject to specific variation in concentration under different pathological conditions<sup>1</sup>.

Since the introduction of moving boundary electrophoresis by Tiselius<sup>2</sup>, and the subsequent use of zone electrophoresis, serum proteins have been fractionated on the basis of their charge at a particular pH. The SAS-I Serum protein kit separates serum proteins into 5 main classes (albumin, alpha 1-globulin, alpha 2-globulin, beta-globulin, and gamma globulin) according to charge in an agarose gel. The proteins are then stained to allow visualisation and quantitative interpretation. Each of the classical electrophoretic zones, with the exception of albumin, normally contains 2 or more components. The relative proportions of these fractions have proven to be useful aids in the diagnosis and prognosis of certain disease states<sup>3,5</sup>.

**WARNINGS AND PRECAUTIONS**

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

**COMPOSITION****1. SAS-I Serum Protein Gel**

Contains agarose in a Tris / Barbitol buffer with thiomersal and sodium azide as preservative. The gel is ready for use as packaged.

**2. Acid Blue Stain Concentrate**

Contains concentrated Acid Blue stain. Dilute the contents of the bottle to 700ml with purified water. Stir overnight and filter before use. Store in a tightly stoppered bottle.

**3. Destain Solution Concentrate**

Dilute the contents of Destain A to 1 litre with purified water. Then add the contents of Destain B and add a further 1 litre of purified water, slowly.

**4. Other Kit Components**

Each kit contains Instructions For Use and sufficient Blotter C to complete 10 gels.

**STORAGE AND SHELF-LIFE****1. SAS-I Serum Protein Gel**

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

**2. Acid Blue Stain**

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain solution is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration of the stain solution.

### 3. Destain Solution

The destain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted destain solution is stable for 6 months at 15...30°C.

#### ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 210200 Sample Applicator Blades (1 x 10)

Cat. No. 210300 Sample Applicator Blades (5 x 10)

Cat. No. 210100 Disposable sample cups (100)

Cat. No. 3100 REP Prep

Drying oven with forced air capable of 60...70°C

Purified water

#### SAMPLE COLLECTION AND PREPARATION

Freshly collected serum is the specimen of choice. Samples can be stored at 15...30°C for up to 4 days, 2...6°C for up to 2 weeks or 6 months at -20°C<sup>6</sup>. Urine and CSF can also be used following a suitable concentration step (50 - 100X). The use of plasma will result in a fibrinogen band between the beta and gamma fractions.

**Interfering Factors:** 1) Haemolysis may cause false elevation in the alpha 2 and beta fractions.  
2) Inaccurate results may be obtained on specimens left uncovered, due to evaporation.

#### STEP-BY-STEP PROCEDURE

1. Pipette 35µl of the sample into the appropriate well of the SAS-I sample tray or disposable sample cups.
  - i) **SAS-I & SAS-I Plus users:** Carefully place the sample tray onto the applicator drawer. Ensure that the tray is pushed firmly down into position.
  - ii) **SAS-3 users:** Carefully locate the sample tray using the sample base locating pins. Ensure that the tray is positioned securely.
2. Remove the gel from the packaging, discard the overlay and:
  - i) **SAS-I users:** place the gel in the SAS-I, agarose side up, aligning the positive and negative sides with the corresponding electrode posts.
  - ii) **SAS-I Plus users:** dispense 400µL of REP Prep onto the heat sink. Place the gel onto the heat sink, agarose side up, aligning the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
  - iii) **SAS-3 users:** place the alignment guide onto the pins and dispense 400µL of REP Prep onto the centre of the chamber. Place the gel into the chamber agarose side up, using the guide, align the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
3. Blot the surface of the gel with a blotter C, discard the blotter.
4.
  - i) **SAS-I users:** attach the electrodes onto the top side of the electrode posts so that they are in contact with the gel blocks.
  - ii) **SAS-I Plus users:** (as above). Place the cover over the gel and electrodes and press firmly for 5 seconds to ensure contact.
  - iii) **SAS-3 users:** attach the electrodes onto the the electrode posts so that they are in contact with the gel blocks.

5. Place two applicator blade assemblies in position on the instrument, (**SAS-3 users:** slot A and 10).
6. Perform the Serum Protein electrophoresis:
  - i) **SAS-1 users:** 80 volts, 20 mins, one application.
  - ii) **SAS-1 Plus users:** Electrophoresis: 100 volts, 18 mins, 20°C, one application
  - iii) **SAS-3 users:**

Step	Time (mm:ss)	Temperature (°C)	Voltage	Other
Load Sample	00:30	21		Speed 1
Apply Sample	00:30	21		Speed 1*
Electrophoresis	18:00	20	100	
Dry	08:00	54		

- \* Use Location 2

**NOTE:** Remove gel blocks prior to drying.

7. Following electrophoresis:
  - i) **SAS-1 Plus users:** remove the cover,
  - ii) **SAS-1 and SAS-1 Plus users:** remove the electrodes and remove both gel blocks using the Gel Block Remover.
8. Attach the gel to the staining chamber holder.
9. Select the Serum Protein test program on the staining unit and following the prompts, Fix, Stain, Destain and Dry the gel.
  - a) **SAS-2 (Auto-Stainer)**

Step	Solution	Time (mm:ss)	Port	Temperature (°C)
Stain	Acid Blue Stain	10:00	6	
Wash	Purified water	01:00	1	
Dry	—	15:00		65
Destain	Destain solution	03:00	2	
Destain	Destain solution	03:00	2	
Wash	Purified water	01:00	1	
Dry	—	10:00		65

- b) **SAS-4 (Auto-Stainer)**

Step	Time (mm:ss)	Temperature (°C)	Other
Stain	04:00		Recirculate ON
Destain	02:00		Recirculate ON
Destain	02:00		Recirculate ON
Dry	12:00	63	

- c) **Manual**

Follow the sequence listed for the SAS-2 Auto-Stainer, using a staining bath for the Stain, Destain and Wash steps, and a Drying Oven with forced air at 60...70°C for the Dry steps.

10. At the end of the staining cycle, remove the gel from the staining chamber. The gel is now ready for examination.

## QUALITY CONTROL

Kemtrol Serum Controls (Cat. No. 7024 and 7025) can be used to verify all phases of the procedure and should be used on each plate run. Refer to the package insert provided for appropriate assay values.

## INTERPRETATION OF RESULTS

It is recommended that any evaluation of the gels is performed against normal values produced for this method in each individual laboratory.

For a complete review of serum protein evaluation, see Ritzmann, S.E, 1982<sup>5</sup>. Studies show that the values are the same for both males and non-pregnant females. Some differences are seen in pregnant females at term and women on oral contraceptives.

Age has some effect on normal levels. Cord blood has a decreased total protein, albumin, alpha 2 and beta fractions; slightly increased alpha 1 and normal or increased gamma fraction (largely of maternal origin). The gamma globulins drop rapidly until about 3 months of age, while other fractions have reached adult levels by this time. Adult levels of the gamma globulins are not reached until 10-16 years of age. The albumin decreases and beta globulin increases over the age of 40.

### 1. Qualitative Evaluation:

The gels may be visually inspected for the presence or absence of particular bands of interest.

### 2. Quantitative Evaluation:

Scan the gels gel, side down, at 595nm.

In either case, an elevation or decrease in particular serum components or the detection of unusual serum components require further investigation. The completed SAS-1 Serum Protein gel is stable for an indefinite period of time.

## LIMITATIONS

Since all electrophoresis procedures are non-linear, it is important to follow these Instructions For Use closely to ensure optimal resolution and reproducible results. Failure to follow these Instructions For Use may affect the results obtained.

## REFERENCE VALUES

Using 20 normal specimens from male and female donors with an age range of 20-59 years, the following normal ranges were obtained (these are presented as a guideline only):

Protein Fraction	Mean (%)	$\pm$ 2SD Range
Albumin	60.4	55.8 - 65.0
Alpha1	3.4	2.2 - 4.6
Alpha2	10.3	8.2 - 12.5
Beta	10.7	7.2 - 14.2
Gamma	15.0	11.5 - 18.6

**PERFORMANCE CHARACTERISTICS****Reproducibility**

	Within Gel (n=16)		Between Gel (n=95)	
	Mean (%)	CV (%)	Mean (%)	CV (%)
Albumin	60.1	1.7	60.4	2.3
Alpha 1	4.2	7.5	4.3	7.3
Alpha 2	11.4	4.1	10.3	7.8
Beta	10.5	2.1	11.2	6.8
Gamma	13.8	3.0	13.8	5.5

**Sensitivity**

The method is sensitive to 0.3 g/L per band, determined as the lowest concentration of protein which was evident as a discrete band on the completed gel.

**Linearity**

The linearity of the method is a function of densitometer specification as well as gel performance. It is recommended that each customer determine the linearity of the method based upon the densitometer in use in the laboratory.

**BIBLIOGRAPHY**

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
2. Tiselius, A. 'A New Apparatus for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Ritzmann, S.E. and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays' in Laboratory Medicine, Harper and Row, Inc., Hagerstown, 1979.
4. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, pages 579-582, 1986.
5. Ritzmann, S.E. (Ed.), Protein Abnormalities Volume I : Physiology of Immunoglobulins - Diagnostic and Clinical Aspects', Allen R. Liss, Inc., New York, 1982.
6. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, (3rd Edition), W.B. Saunders Co., Philadelphia, page 524, 1995.

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**LITERATUR**

- Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
- Tiselius, A. 'A New Approach for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
- Ritzmann, S.E. and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays' in Laboratory Medicine, Harper and Row, Inc., Hagerstown, 1979.
- Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, Seite 579-582, 1986.
- Ritzmann, S.E. (Ed.), Protein Abnormalities Volume I : 'Physiology of Immunoglobulins - Diagnostic and Clinical Aspects', Allen R. Liss, Inc., New York, 1982.
- Tietz, N.W. (Ed.), Textbook of Clinical Chemistry (3rd Edition), W.B. Saunders Co., Philadelphia, Seite 524, 1995.

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**BIBLIOGRAFIA**

- Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
- Tiselius, A. 'A New Apparatus for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
- Ritzmann, S.E. and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays' in Laboratory Medicine, Harper and Row, Inc., Hagerstown, 1979.
- Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, pages 579-582, 1986.
- Ritzmann, S.E. (Ed.), Protein Abnormalities Volume I : Physiology of Immunoglobulins - Diagnostic and Clinical Aspects', Allen R. Liss, Inc., New York, 1982.
- Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, (3rd Edition), W.B. Saunders Co., Philadelphia, page 524, 1995.

# helena | BioSciences Europe

[www.helena-biosciences.com](http://www.helena-biosciences.com)

---

Helena BioSciences Europe  
Queensway South  
Team Valley Trading Estate  
Gateshead  
Tyne and Wear  
NE11 0SD

tel: +44 (0) 191 482 8440  
fax: +44 (0) 191 482 8442

email: [info@helena-biosciences.com](mailto:info@helena-biosciences.com)