## **SAS-I Acid Hb-12**

### **Instructions For Use REF 201000**

**SAS-I Hb-12 Acide**  
Fiche technique

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**Hb-12 Ácida SAS-I**  
Instrucciones de uso

### **Contents**

English .....	1
Français .....	8
Deutsch .....	15
Italiano .....	22
Español .....	29

**INTENDED PURPOSE**

The SAS-I Acid Hb-12 kit is intended for the separation of human haemoglobins by agarose gel electrophoresis.

Haemoglobins (Hb) are a group of proteins whose chief functions are to transport oxygen from the lungs to the tissues and carbon dioxide in the reverse direction. They are composed of polypeptide chains called globin, and iron protoporphyrin haem groups. A specific sequence of amino acids constitutes each of the four polypeptide chains. Each normal haemoglobin molecule contains one pair of alpha and one pair of non-alpha chains. In normal adult haemoglobin (HbA), the non-alpha chains are called beta. The non-alpha chains of foetal haemoglobin are called gamma. A minor (3%) haemoglobin fraction called HbA<sub>2</sub> contains alpha and delta chains. Two other chains are formed in the embryo.

The major haemoglobin in the erythrocytes of the normal adult is HbA and there are small amounts of HbA<sub>2</sub> and HbF. In addition, over 400 mutant haemoglobins are now known, some of which may cause serious clinical effects, especially in the homozygous state or in combination with another abnormal haemoglobin. Wintrobe<sup>1</sup> divides the abnormalities of haemoglobin synthesis into three groups.

1. Production of an abnormal protein molecule (e.g. sickle cell anaemia).
2. Reduction in the amount of normal protein synthesis (e.g. thalassaemia).
3. Developmental anomalies (e.g. hereditary persistence of foetal haemoglobin (HPFH)).

The two mutant haemoglobins most commonly seen are HbS and HbC. Hb Lepore, HbE, HbG-Philadelphia, HbD-Los Angeles, and HbO-Arab may be seen less frequently<sup>2</sup>.

Electrophoresis is generally considered the best method for separating and identifying haemoglobinopathies. The protocol for haemoglobin electrophoresis involves step-wise use of two systems<sup>3,4</sup>. Initial electrophoresis is performed in alkaline buffers. However, because of the electrophoretic similarity of many structurally different haemoglobins, the evaluation must be supplemented by acid buffer electrophoresis which measures a property other than electrical charge.

This method is based on the complex interactions of the haemoglobin with an acid electrophoretic buffer and the agarose support. The SAS-I Acid Hb-12 procedure is a simple procedure requiring minute quantities of haemolysate to provide complementary evidence (along with the results from SAS-I Alkaline Hb-12 analysis) of the presence of HbS, HbC and HbF as well as several other abnormal haemoglobins.

**WARNINGS AND PRECAUTIONS**

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

**COMPOSITION****1. SAS-I Acid Hb Gel**

Contains agarose in a Citrate / Maleate buffer with sodium azide as preservative. The gel is ready for use as packaged.

## 2. Acid Violet Stain Concentrate

Contains concentrated Acid Violet stain. Dilute the contents of the bottle to 700ml with purified water. Stir overnight and filter before use. Store in a tightly stoppered bottle.

## 3. Haemoglobin Lysing Agent

Contains Triton X-100 in purified water with potassium cyanide and thiomersal as preservative. The Lysing Agent is ready for use as packaged.

## 4. Other Kit Components

Each kit contains Instructions For Use and sufficient Blotter C and Soak-away Blotters to complete 10 gels.

## STORAGE AND SHELF-LIFE

### 1. SAS-I Acid Hb Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

### 2. Acid Violet Stain

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain solution is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration.

### 3. Haemoglobin Lysing Agent

The Lysing Agent should be stored at 15...30°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

## ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 210200 Sample Applicator Blades (1 x 10)

Cat. No. 210300 Sample Applicator Blades (5 x 10)

Cat. No. 210100 Disposable sample cups (100)

Cat No. 310907 SPI00 Sample Tray

Cat. No. 3100 REP Prep

Cat. No. 5331 AFSC Hemo Control

Cat. No. 5330 AFSA<sub>2</sub> Hemo Control

Cat. No. 5329 ASA<sub>2</sub> Hemo Control

Cat. No. 5328 AA<sub>2</sub> Hemo Control

Drying oven with forced air capable of 60...70°C

Fixative / Destain Solution: Mix 50ml of glacial acetic acid and 950ml of purified water. Store in a tightly stoppered bottle.

Saline solution (0.85% NaCl)

Purified water

## SAMPLE COLLECTION AND PREPARATION

Freshly collected EDTA or heparin anticoagulated blood is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 1 week. For optimal results, saline washed red cells should be used to prepare lysates. This removes possible interference from plasma proteins.

a) Mix 200µL of well mixed whole blood with 1000µL of saline solution.

b) Centrifuge to sediment the red cells.

- c) Remove 1000 $\mu$ L of the supernatant and discard.
- d) Add a further 1000 $\mu$ L of saline solution and mix well.
- e) Repeat steps b-d x2.
- f) Following the final centrifugation, remove 1000 $\mu$ L of supernatant and treat the remaining sample as whole blood, or remove all of the supernatant and treat the remaining sample as washed packed cells.

Dilute each patient sample / control to a haemoglobin concentration of 1.0-2.0g/dL with Haemoglobin Lysing Agent.

### **STEP-BY-STEP PROCEDURE**

1. Pipette 35 $\mu$ l of the sample into the appropriate well of the sample tray or disposable sample cups.
  - i) **SAS-I & SAS-I Plus users:** Use SAS-I sample tray. Carefully place the sample tray onto the applicator drawer. Ensure that the tray is pushed firmly down into position.
  - ii) **SAS-3 users:** Use SP100 sample tray. Carefully locate the sample tray using the sample base locating pins. Ensure that the tray is positioned securely.
2. Remove the gel from the packaging and:
  - i) **SAS-I users:** place the gel in the SAS-I, agarose side up, aligning the positive and negative sides with the corresponding electrode posts.
  - ii) **SAS-I Plus users:** dispense 400 $\mu$ L of REP Prep onto the heat sink. Place the gel onto the heat sink, agarose side up, aligning the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
  - iii) **SAS-3 users:** place the alignment guide onto the pins and dispense 400 $\mu$ L of REP Prep onto the centre of the chamber. Place the gel into the chamber agarose side up, using the guide, align the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel. Position one soak-away blotter on the outside of the cathode electrode and one soak-away blotter on the outside of the anode electrode to act as excess fluid blotters. Ensure the soak-aways are in contact with the gel.
3. Blot the surface of the gel with a blotter C, discard the blotter.
4.
  - i) **SAS-I users:** attach the electrodes onto the top side of the electrode posts so that they are in contact with the gel blocks.
  - ii) **SAS-I Plus users:** (as above). Place the cover over the gel and electrodes and press firmly for 5 seconds to ensure contact.
  - iii) **SAS-3 users:** attach the electrodes onto the electrode posts so that they are in contact with the gel blocks.
5. Place I applicator blade assembly into the lower position on the instrument, (**SAS-3 users:** Slot 9).
6. Perform the Acid Haemoglobin electrophoresis:
  - i) **SAS-I users:** 60 volts, 22 mins, 2 applications.
  - ii) **SAS-I Plus users:** Electrophoresis: 60 volts, 22 mins, 18°C, 2 applications.

iii) **SAS-3 users:**

Step	Time (mm:ss)	Temperature (°C)	Voltage	Other
Load Sample	00:30	18		Speed 6
Apply Sample	00:30	18		Speed 6*
<b>Repeat steps 'load sample, apply sample' for a total of 2 applications.</b>				
Electrophoresis	14:00	18	100	
Dry	08:00	62		

\* Use Location 2

**NOTE:** Remove gel blocks prior to drying.

7. Following electrophoresis:

i) **SAS-1 Plus users:** remove the cover.

ii) **SAS-3 users:** remove the soak away blotters.

iii) **SAS-1 and SAS-1 Plus users:** remove both gel blocks using the Gel Block Remover.

8. Attach the gel to the staining chamber holder.

9. Select the Acid Haemoglobin test program on the staining unit and, following the prompts, Fix, Stain, Destain and Dry the gel.

a) **SAS-2 (Auto-Stainer)**

Step	Solution	Time (mm:ss)	Port	Temperature (°C)
Fix	Fixative / Destain solution	05:00	4	
Stain	Acid Violet	15:00	5	
Destain	Fixative / Destain solution	00:10	4	
Dry	—	15:00		65
Destain	Fixative / Destain solution	02:00	4	
Destain	Fixative / Destain solution	02:00	4	
Wash	Purified water	01:00	1	
Dry	—	15:00		65

b) **SAS-3 users (Manual)**

Follow the sequence listed for the SAS-2 Auto-Stainer, using a staining bath for the Fix, Stain, Destain and Wash steps, and a Drying Oven with forced air at 60...70°C for the Dry steps.

10. At the end of the staining cycle, remove the gel from the staining chamber. The gel is now ready for examination.

## INTERPRETATION OF RESULTS

Qualitative Evaluation : The possible identity of the haemoglobin types present in the samples can be determined by visual evaluation of the completed gel. The Hemo controls provide a marker for band identification.

Figure 1 shows the possible identity of the most commonly encountered haemoglobin types.


FIG. 1

Most haemoglobin variants cause no discernible clinical symptoms, so are of interest primarily to research scientists. Variants are clinically important when their presence leads to sickling disorders, thalassaemia syndromes, life long cyanosis, haemolytic anaemias or erythrocytosis, or if the heterozygote is of sufficient prevalence to warrant genetic counselling. The combinations of HbS-S, HbS-D-Los Angeles, and HbS-O Arab lead to serious sickling disorders<sup>2</sup>. Several variants including HbH, E-Fort Worth and Lepore cause a thalassaemic blood picture<sup>2</sup>.

The two variant hemoglobins of greatest importance in terms of frequency and pathology are HbS and HbC<sup>2</sup>. Sickle cell anaemia (HbSS) is a cruel and lethal disease. It first manifests itself at about 5-6 months of age. The clinical course presents agonising episodes of pain and temperature elevations with anaemia, listlessness, lethargy, and infarct in virtually all organs of the body. The individual with homozygous HbCC suffers mild haemolytic anaemia which is attributed to the precipitation or crystallization of HbC within the erythrocytes. Cases of HbSC disease are characterised by haemolytic anaemia that is milder than sickle-cell anaemia.

The thalassaemias are a group of haemoglobin disorders characterised by hypochromia and microcytosis due to the diminished synthesis of one globin chain (the a or b) while synthesis of the other chain proceeds normally<sup>9,10</sup>. This unbalanced synthesis results in unstable globin chains. These precipitate within the red cell, forming inclusion bodies that shorten the life span of the cell. In a-thalassaemias, the a chains are diminished or absent, and in the b-thalassaemia, the b chains are affected. Another quantitative disorder of haemoglobin synthesis, hereditary persistent foetal haemoglobin (HPFH), represents a genetic failure of the mechanisms that turn off gamma chain synthesis at about four months after birth, which results in a continued high percentage of HbF. It is a more benign condition than the true thalassaemias and persons homozygous for HPFH have normal development, are asymptomatic and have no anemia<sup>10</sup>.

The most common haemoglobin abnormalities:

#### **Sickle Cell Trait**

This is a heterozygous state showing HbA and HbS and a normal amount of HbA<sub>2</sub> on cellulose acetate. Results on citrate agar show haemoglobins in the HbA and HbS migratory positions (zones).

#### **Sickle Cell Anaemia**

This is a homozygous state showing almost exclusively HbS, although a small amount of HbF may also be present.

#### **Sickle-C Disease**

This is a heterozygous state demonstrating HbS and HbC.

### **Sickle Cell-Thalassaemia Disease**

This condition shows HbA, HbF, HbS, and HbA<sub>2</sub>.

In Sickle Cell  $\beta^0$ -Thalassaemia HbA is absent.

In Sickle Cell  $\beta^+$ -Thalassaemia HbA is present in reduced quantities.

### **Thalassaemia-C Disease**

This condition shows HbA, HbF and HbC.

### **C Disease**

This is a homozygous state showing almost exclusively HbC.

**Thalassaemia Major** This condition shows HbF, HbA and HbA<sub>2</sub>.

## **LIMITATIONS**

Some abnormal haemoglobins have similar electrophoretic mobilities and must be differentiated by other methodologies.

Further testing required:

1. Globin chain analysis (both acid and alkaline) and structural studies may be necessary in order to positively identify some of the more rare haemoglobins.
2. Anion exchange column chromatography is the most accurate method for quantitating HbA<sub>2</sub>. Helena BioSciences Sickle-Thal Quik Column Method (Cat. No. 5334) for quantitation of HbA<sub>2</sub> in the presence of HbS, or the Helena BioSciences Beta-Thal HbA<sub>2</sub> Quik Column Procedure (Cat. No. 5341) are recommended. HbA<sub>2</sub> quantitation is one of the most important diagnostic tests in the diagnosis of B-thalassaemia trait.
3. Low levels of HbF (1-10%) may be accurately quantitated by radial immunodiffusion using the Helena BioSciences HbF-QuikPlate Procedure (Cat. No. 9325).

## **REFERENCE VALUES**

At birth, the majority of haemoglobin in the erythrocytes of the normal individual is foetal haemoglobin, HbF. Some of the major adult haemoglobin, HbA, and a small amount of HbA<sub>2</sub>, are also present. At the end of the first year of life and through adulthood, the major haemoglobin present is HbA with up to 3.7% HbA<sub>2</sub> and less than 2% HbF.

## **PERFORMANCE CHARACTERISTICS**

### **a) Reproducibility**

The SAS-I Acid Hb gel is a qualitative system for the identification of haemoglobin bands. Using a control material containing Hb's A, S and A<sub>2</sub>, the same band patterns were seen within a single gel and between different gels. No bands were missing and no additional bands were observed between applications.

### **b) Sensitivity**

0.08 g/dL per band, determined as the lowest concentration of haemoglobin which was evident as a discrete band on the completed gel.

### **c) Linearity**

4 g/dL per band, based upon the maximum concentration of haemoglobin per band which allows separated bands to be satisfactorily resolved from each other without protein overloading. This procedure is not intended for densitometric scanning.

**BIBLIOGRAPHY**

1. Wintrobe, Maxwell M., Clinical Hematology, 6th Edition, Lea and Febiger, Philadelphia, 1967. pages 145-167.
2. Fairbanks, V.F., 'The Nomenclature and Taxonomy of Hemoglobin Variants', Diagnostic Medicine, Nov/Dec., 53-58, 1980.
3. Schneider, R.G., Hightower, B.J. and Barwick, R.C., 'Laboratory Identification of the Hemoglobins', Lab Management, 1981; August: 29-43.
4. Center for Disease Control, Laboratory Methods for Detecting Hemoglobinopathies, U.S. Department of Health and Human Services/Public Health Service, 1984.
5. Schneider, R.G., 'Methods for Detection of Hemoglobin Variants and Hemoglobinopathies in the Routine Clinical Laboratory', CRC Critical Reviews in Clinical Laboratory Sciences, 1978.
6. Schneider, R.G., Hightower, B., Hosty, T.S., Ryder, H., Tomlin, G., Atkins, R., Brimhall, B., and Jones, R.T., 'Abnormal Hemoglobins in a Quarter Million People', Blood, 1976; 48(5) : 629-637.
7. Huisman, T.H.J. and Schroeder, W.A., 'New Aspects of the Structure, Function and Synthesis of Hemoglobins'. CRC Press, Cleveland, 1971.
8. Schmidt, R.M., Huisman, T.H.J., and Lehmann, H., The Detection of Hemoglobinopathies. CRC Press, Cleveland, 1974.
9. Weatherall, D.J. and Clegg, J.B., The Thalassemia Syndromes, Blackwell Scientific Publications, Oxford, 1972.
10. Lehman, H. and Huntsman, R.G., Man's Haemoglobins, J.B. Lippincott Co., Philadelphia, 1974.

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**LITERATUR**

1. Wintrobe, Maxwell M., Clinical Hematology, 6th Edition, Lea and Febiger, Philadelphia, 1967. Seite 145-167
2. Fairbanks, V.F., 'The Nomenclature and Taxonomy of Hemoglobin Variants', Diagnostic Medicine, Nov/Dec., 53-58, 1980.
3. Schneider, R.G., Hightower, B.J. and Barwick, R.C., 'Laboratory Identification of the Hemoglobins', Lab Management, 1981; August: 29-43.
4. Center for Disease Control, Laboratory Methods for Detecting Hemoglobinopathies, U.S. Department of Health and Human Services/Public Health Service, 1984.
5. Schneider, R.G., 'Methods for Detection of Hemoglobin Variants and Hemoglobinopathies in the Routine Clinical Laboratory', CRC Critical Reviews in Clinical Laboratory Sciences, 1978.
6. Schneider, R.G., Hightower, B., Hosty, T.S., Ryder, H., Tomlin, G., Atkins, R., Brimhall, B., and Jones, R.T., 'Abnormal Hemoglobins in a Quarter Million People', Blood, 1976; 48(5) : 629-637.
7. Huisman, T.H.J. and Schroeder, W.A., 'New Aspects of the Structure, Function and Synthesis of Hemoglobins'. CRC Press, Cleveland, 1971.
8. Schmidt, R.M., Huisman, T.H.J., and Lehmann, H., The Detection of Hemoglobinopathies. CRC Press, Cleveland, 1974.
9. Weatherall, D.J. and Clegg, J.B., The Thalassemia Syndromes, Blackwell Scientific Publications, Oxford, 1972.
10. Lehman, H. and Huntsman, R.G., Man's Haemoglobins, J.B. Lippincott Co., Philadelphia, 1974.

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